

This document is provided for product evaluation purposes only. It is not intended to be used in place of the package insert shipped with the product.



<http://www.zeptometrix.com/>

1-800-274-5487



## **NATtrol™ BK Virus Linearity Panel**

ZMC Catalog #: NATBK-Lin

**FOR RESEARCH USE ONLY**  
Not for *in vitro* Diagnostic Use

### **Product Description**

NATtrol™ BK Linearity Panel (NATBK-Lin) is formulated with intact BK virus particles that have been chemically modified to render them non-infectious and refrigerator stable\*. Each linearity panel contains 6 x 0.25 mL vials of each member with concentrations spanning the clinically significant range for BK infected individuals. These controls are supplied in a purified serum protein matrix that mimics the composition of a true clinical specimen. Panel members are quantitated using an in-house real-time PCR assay targeting the BK Virus T antigen gene. Copies/mL were calibrated against commercially available purified BK virus DNA.

\*NATtrol™ Patents Pending.

### **Intended Use**

- The NATtrol™ BK Linearity Panel (NATBK-Lin) is designed to evaluate the recovery, accuracy and sensitivity of nucleic acid tests (NAT) for the determination of BK virus DNA. NATBK-Lin can also be used for assay validation, training of laboratory personnel, QA/QC and educational purposes.
- The NATBK-Lin panel contains intact organisms and should be run in a manner identical to that used for clinical specimens.
- **All panel members must be extracted prior to amplification.**

The table below lists the targeted copies / ml concentrations for each member in the linearity panel.

Panel Member	BK Virus DNA copies/ml
100 M	10 <sup>8</sup>
10 M	10 <sup>7</sup>
1 M	10 <sup>6</sup>
100K	10 <sup>5</sup>
10k	10 <sup>4</sup>
Negative	0

- The purified serum protein matrix was sourced from licensed U.S. blood banks and screened negative for HIV 1&2 Ab, HBsAg, HTLV I&II Ab, HCV Ab, HIV RNA, HBV DNA and HCV RNA using FDA cleared kits at the individual donor level.

***DO NOT USE IN HUMANS OR AS A CLINICAL DIAGNOSTIC.***

These products are intended for research, product development or manufacturing use only. These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

**Precautions**

- Use Universal Precautions when handling NATBK-Lin.
- To avoid cross-contamination, use separate pipette tips for all reagents.



This product was manufactured in a facility whose Quality Management System is certified as being in compliance with ISO 9001:2008 and ISO 13485:2003 standards

ZeptoMetrix Corporation  
872 Main Street  
Buffalo, New York 14202  
Office Phone: 716-882-0920  
Fax: 716-882-0959  
Customer Service: 800-274-5487  
www.zeptometrix.com  
E-mail: support@zeptometrix.com  
custserv@zeptometrix.com