



**HBV Longitudinal Panel  
Part Number: HBV6536**

*These reagents are not a substitute for the mandatory positive and negative control reagents provided with licensed test kits.  
For Research Use Only. Not for use in diagnostic procedures.*

**\*Donor is co-infected with Hepatitis C.**

| PANEL NUMBER | DATE OF DRAW | ORTHO VITROS ECi ANTI-HIV ½ | ORTHO VITROS ECi ANTI-HCV | ABBOTT HBsAg   | ORTHO VITROS ECi HBsAg | ABBOTT ANTI-HBc TOTAL | ORTHO VITROS ECi ANTI-HBc TOTAL | ABBOTT ANTI-HBc IgM |
|--------------|--------------|-----------------------------|---------------------------|----------------|------------------------|-----------------------|---------------------------------|---------------------|
|              |              | TEST DATE: 8/21/2014        | TEST DATE: 8/21/2014      | TEST DATE: UNK | TEST DATE: 8/21/2014   | TEST DATE: UNK        | TEST DATE: 8/21/2014            | TEST DATE: UNK      |
| 6536-01      | 11/13/01     | 0.10                        | 34.9                      | >36.4          | 8340                   | 0.04                  | 0.09                            | >4.7                |
| 6536-02      | 11/26/01     | 0.07                        | 34.8                      | >36.4          | 2060                   | 0.03                  | 0.14                            | >4.7                |
| 6536-03      | 12/11/01     | 0.09                        | 35.2                      | >36.4          | 44.9                   | 0.07                  | 0.11                            | >4.7                |
| 6536-04      | 12/25/01     | 0.08                        | 34.7                      | 5.6            | 3.57                   | 0.05                  | 0.04                            | >4.7                |
| 6536-05      | 01/08/02     | 0.11                        | 36.6                      | 0.4            | 0.13                   | 0.11                  | 0.04                            | 3.9                 |
| 6536-06      | 01/22/02     | 0.11                        | 35.6                      | 0.7            | 0.10                   | 0.09                  | 0.02                            | 3.1                 |
| 6536-07      | 02/05/02     | 0.10                        | 34.8                      | 0.2            | 0.10                   | 0.06                  | 0.03                            | 2.5                 |
| 6536-08      | 02/19/02     | 0.09                        | 34.9                      | 0.2            | 0.09                   | 0.07                  | 0.02                            | 2.4                 |
| 6536-09      | 03/05/02     | 0.10                        | 35.0                      | 0.2            | 0.09                   | 0.07                  | 0.03                            | 2.1                 |
| 6536-10      | 03/19/02     | 0.11                        | 35.4                      | 0.2            | 0.09                   | 0.07                  | 0.02                            | 1.8                 |
| 6536-11      | 04/02/02     | 0.09                        | 35.2                      | 0.1            | 0.09                   | 0.07                  | 0.02                            | 2.0                 |
| 6536-12      | 04/16/02     | 0.10                        | 35.2                      | 0.2            | 0.10                   | 0.06                  | 0.02                            | 1.7                 |

| PANEL NUMBER | ORTHO VITROS ECi ANTI-HBc IgM | DIASORIN ANTI-HBe | ORTHO VITROS ECi ANTI-HBe | ABBOTT ANTI-HBs QUANT | ORTHO VITROS ECi ANTI-HBs QUANT | ORTHO VITROS ECi HBeAg | DIASORIN HBeAg |
|--------------|-------------------------------|-------------------|---------------------------|-----------------------|---------------------------------|------------------------|----------------|
|              | TEST DATE: 8/21/2014          | TEST DATE: UNK    | TEST DATE: 8/21/2014      | TEST DATE: UNK        | TEST DATE: 8/21/2014            | TEST DATE: 8/21/2014   | TEST DATE: UNK |
| 6536-01      | 3.70                          | 0.3               | 0.97                      | 0.3                   | 1.6                             | 0.40                   | 0.31           |
| 6536-02      | 3.66                          | 0.4               | 1.45                      | 0.3                   | 1.7                             | 0.14                   | 0.31           |
| 6536-03      | 3.02                          | 0.1               | 9.54                      | 0.4                   | 2.5                             | 0.14                   | 0.34           |
| 6536-04      | 2.37                          | 0.07              | 26.8                      | 0.3                   | 2.7                             | 0.12                   | 0.26           |
| 6536-05      | 1.66                          | 0.05              | 20.6                      | 0.5                   | 3.5                             | 0.20                   | 0.07           |
| 6536-06      | 1.43                          | 0.05              | 26.7                      | 0.7                   | 5.2                             | 0.21                   | 0.16           |
| 6536-07      | 1.28                          | 0.03              | 27.7                      | 0.9                   | 8.5                             | 0.19                   | 0.09           |
| 6536-08      | 1.21                          | 0.04              | 8.28                      | 1.3                   | 12.1                            | 0.19                   | 0.23           |
| 6536-09      | 1.10                          | 0.04              | 26.9                      | 1.8                   | 19.5                            | 0.19                   | 0.16           |
| 6536-10      | 1.01                          | 0.05              | 25.6                      | 2.2                   | 20.3                            | 0.17                   | 0.26           |
| 6536-11      | 1.04                          | 0.05              | 25.7                      | 2.2                   | 33.1                            | 0.18                   | 0.16           |
| 6536-12      | 0.99                          | 0.05              | 19.2                      | 2.1                   | 26.6                            | 0.19                   | 0.24           |

**Notes on Data:**

1. This panel is for **Research Use Only. Not for use in diagnostic procedures.**
2. Data is generated at ZeptoMetrix™ LLC and by our Laboratory Partners in the United States and Europe.
3. All specimens collected in this longitudinal series are unadulterated 4% sodium citrated plasma samples collected from a single donor in the United States.
4. No preservatives have been added.

ZeptoMetrix LLC • 25 Kenwood Circle, Franklin, MA 02038 • Tel (508) 553-5800 • Fax (508) 520-1525  
This Product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.  
For Customer Support, please visit [www.zeptometrix.com](http://www.zeptometrix.com) or email [zepto.customerservice@antylia.com](mailto:zepto.customerservice@antylia.com)

PIHBV6536 Rev02  
Effective Date: 09/07/2021