



HBV VACCINE PANEL

(Catalog No. 3893)

For Investigational Research Use Only

Application Notes: This 10 member panel was collected from a 29 year old Caucasian Female.

PANEL NUMBER	DATE OF DRAW	HBV RESULT
3893-01	12/3/96	0.2125
3893-02	1/17/97	0.175
3893-03	1/31/97	0.1875
3893-04	2/14/97	0.2375
3893-05	2/28/87	0.5625
3893-06	3/26/97	0.175
3893-07	4/24/97	0.256
3893-08	5/29/97	0.632
3893-09	6/26/97	0.545
3893-10	7/31/97	6.039