



HBV VACCINE PANEL

(Catalog No. 3894)

For Investigational Research Use Only

Application Notes: This 10 member panel was collected from a 21 year old Caucasian Female.

| PANEL NUMBER | DATE OF DRAW | HBV RESULT |
|--------------|--------------|---------------|
| 3894-01 | 12/3/96 | 0.25 |
| 3894-02 | 1/17/97 | 0.25 |
| 3894-03 | 1/31/97 | 2.725 |
| 3894-04 | 2/14/97 | 1.4125 |
| 3894-05 | 2/28/97 | 3.8625 |
| 3894-06 | 3/26/97 | 15.85 |
| 3894-07 | 4/24/97 | 13.463 |
| 3894-08 | 5/29/97 | 11.71 |
| 3894-09 | 6/26/97 | 7.792 |
| 3894-10 | 7/31/97 | 21.565 |