



HBV LONGITUDINAL PANEL

(Catalog No. 6506)

For Investigational Research Use Only

Application Notes: This 14 member panel was collected from a 32 year old Caucasian Male.

PANEL NUMBER	DATE OF DRAW	DATE CRF COMPLETED	DATE OF BIRTH	DODX	SIGNS AND SYMPTOMS	RISK GROUP	MODE OF TRANSMISSION	ALT	ALT DATE	BILIRUBIN	BILIRUBIN DATE	SGOT	SGOT DATE
6506-01	3/30/01	4/3/01	8/20/69	3/30/01	JAUNDICE/ANOREXIA/ FATIGUE/HEPATOMEGALY/DARK URINE/MALAISE	NO KNOWN RISK	NOT DETERMINED	570	3/30/01	134	3/30/01	310	3/30/01
6506-02	4/12/01	4/13/01	8/20/69	3/30/01	JAUNDICE/FATIGUE/HEPATOMEGALY/DARK URINE/MALAISE	NO KNOWN RISK	NOT DETERMINED	269	4/12/01	41.9	4/12/01	151	4/12/01
6506-03	4/26/01	4/27/01	8/20/69	3/30/01	JAUNDICE/HEPATOMEGALY/ DARK URINE	NO KNOWN RISK	NOT DETERMINED	36	4/26/01	75	4/26/01	20.9	4/26/01
6506-04	5/10/01	5/11/01	8/20/69	3/30/01	HEPATOMEGALY/ DARK URINE	NO KNOWN RISK	NOT DETERMINED	42	5/10/01	60	5/10/01	34	5/10/01
6506-05	5/24/01	5/25/01	8/20/69	3/30/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	38	5/24/01	24	5/24/01	36	5/24/01
6506-06	6/7/01	6/8/01	8/20/69	3/30/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	31.5	6/7/01	26.1	6/7/01	33.7	6/7/01
6506-07	6/21/01	6/25/01	8/20/69	3/30/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	54	6/22/01	14.5	6/22/01	32	6/22/01
6506-08	7/5/01	7/11/01	8/20/69	3/30/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	12	7/6/01	10.6	7/6/01	22	7/6/01
6506-09	7/19/01	7/27/01	8/20/69	3/30/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	6.6	7/20/01	11.1	7/20/01	35.3	7/20/01
6506-10	8/2/01	8/6/01	8/20/69	3/30/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	18.2	8/3/01	27.5	8/3/01	42.3	8/3/01
6506-11	8/18/01	8/21/01	8/20/69	3/30/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	21.0	8/19/01	13.2	8/19/01	23.5	8/19/01
6506-12	8/30/01	9/1/01	8/20/69	3/30/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	14.7	8/31/01	14.8	8/31/01	31.2	8/31/01
6506-13	11/27/01	11/29/01	8/20/69	3/30/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	13.6	11/28/01	9.7	11/28/01	26.2	11/28/01
6506-14	12/25/01	12/27/01	8/20/69	3/30/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	8	12/26/01	7.7	12/26/01	16	12/26/01



HBV LONGITUDINAL PANEL

(Catalog No. 6506)

For Investigational Research Use Only

Application Notes: This 14 member panel was collected from a 32 year old Caucasian Male.

PANEL NUMBER	MEDICATIONS	CONCURRENT MEDICAL CONDITIONS	ABBOTT HBSAG S/CO	ABBOTT HBC TOTAL S/CO	ABBOTT HBC IGM S/CO	DIASORIN HBEAG S/CO	DIASORIN ANTI-HBE S/CO	ABBOTT ANTI-HBS S/CO	ABBOTT ANTI-HIV 1/2 S/CO	ABBOTT HCV S/CO
6506-01	AGUA MINERALE/JUICEMULTIVITAMIN/SORBENT/SOL. ROBOXINI/VITAMIN AEVIT/ACIDIASCORBINIC	NONE	>40.8	0.1	>6.9	0.2	0.08	0.4	0.2	0.2
6506-02	AGUA MINERALE/JUICE/SORBENT/VITAMIN AEVIT/VITAMIN ACIDI ASCORBINICI	NONE	>40.8	0.05	>6.9	0.2	0.03	0.4	N/A	N/A
6506-03	HEPATIC TEA/ DUOVIT	NONE	5.7	0.08	>6.9	0.1	0.02	0.4	N/A	N/A
6506-04	NONE	NONE	0.02	0.30	>6.9	0.1	0.03	0.5	N/A	N/A
6506-05	NONE	NONE	0.04	0.08	>6.9	0.3	0.02	0.5	N/A	N/A
6506-06	NONE	NONE	0.02	0.10	6.1	0.1	0.03	0.8	N/A	N/A
6506-07	NONE	NONE	0.02	0.08	5.2	0.2	0.01	1.1	N/A	N/A
6506-08	NONE	NONE	0.02	0.09	4.7	0.2	0.02	2.0	N/A	N/A
6506-09	DIET N5	NONE	0.08	0.04	1.14	0.01	0.3	1.15	N/A	N/A
6506-10	DIET N5	NONE	0.02	0.03	3.9	0.01	0.02	6.3	N/A	N/A
6506-11	DIET N5	NONE	0.13	0.03	3.6	0.03	0.03	6.5	N/A	N/A
6506-12	DIET N5	NONE	0.08	0.02	3.3	0.01	0.03	6.6	N/A	N/A
6506-13	DIET N5	NONE	0.05	0.07	1.74	0.01	0.06	5.2	N/A	N/A
6506-14	DIET N5	NONE	0.04	0.06	2.1	0.03	0.04	5.7	N/A	N/A