



# HBV LONGITUDINAL PANEL

(Catalog No. 6507)

For Investigational Research Use Only

**Application Notes:** This 14 member panel was collected from a 17 year old Caucasian Male.

PANEL NUMBER	DATE OF DRAW	DATE CRF COMPLETED	DATE OF BIRTH	DODX	SIGNS AND SYMPTOMS	RISK GROUP	MODE OF TRANSMISSION	ALT	ALT DATE	BILIRUBIN	BILIRUBIN DATE	SGOT	SGOT DATE
6507-01	4/16/01	4/18/01	8/5/84	4/16/01	JAUNDICE/ANOREXIA/FATIGUE/HEPATOMEGALY/NAUSEA/VOMITING/DARK URINE/MALISE/ABDOMINAL PAIN	IVDA	TRANSFUSION	1620	4/16/01	101	4/16/01	846	4/16/01
6507-02	4/30/01	5/3/01	8/5/84	4/16/01	ASYMPTOMATIC	IVDA	TRANSFUSION	394	4/30/01	18.9	4/30/01	241	4/30/01
6507-03	5/14/01	5/21/01	8/5/84	4/16/01	JAUNDICE/FATIGUE/HEPATOMEGALY	IVDA	TRANSFUSION	290	5/15/01	30.0	5/15/01	173.0	5/15/01
6507-04	5/28/01	5/30/01	8/5/84	4/16/01	HEPATOMEGALY	IVDA	TRANSFUSION	268	5/29/01	29.0	5/29/01	154	5/29/01
6507-05	6/11/01	6/12/01	8/5/84	4/16/01	ASYMPTOMATIC	IVDA	TRANSFUSION	26.5	6/11/01	7.9	6/11/01	37.7	6/11/01
6507-06	6/25/01	7/3/01	8/5/84	4/16/01	ASYMPTOMATIC	IVDA	TRANSFUSION	34.8	6/26/01	16.7	6/26/01	42.6	6/26/01
6507-07	7/9/01	7/13/01	8/5/84	4/16/01	ASYMPTOMATIC	IVDA	TRANSFUSION	89.2	7/10/01	15.1	7/10/01	51.8	7/10/01
6507-08	7/23/01	7/27/01	8/5/84	4/16/01	ASYMPTOMATIC	IVDA	TRANSFUSION	49	7/24/01	9.5	7/24/01	32	7/24/01
6507-09	8/6/01	8/16/01	8/5/84	4/16/01	ASYMPTOMATIC	IVDA	TRANSFUSION	16.6	8/15/01	10.1	8/15/01	11.3	8/15/01
6507-10	8/21/01	8/26/01	8/5/84	4/16/01	ASYMPTOMATIC	IVDA	TRANSFUSION	37.5	8/23/01	14.2	8/23/01	38.1	8/23/01
6507-11	9/3/01	9/7/01	8/5/84	4/16/01	ASYMPTOMATIC	IVDA	TRANSFUSION	143	9/7/01	16.2	9/7/01	71.4	9/7/01
6507-12	9/17/01	9/21/01	8/5/84	4/16/01	ASYMPTOMATIC	IVDA	TRANSFUSION	40.0	9/20/01	11.0	9/20/01	39.6	9/20/01
6507-13	11/27/01	11/29/01	8/5/84	4/16/01	ASYMPTOMATIC	IVDA	TRANSFUSION	51.2	11/28/01	10.0	11/28/01	45.9	11/28/01
6507-14	12/11/01	12/14/01	8/5/84	4/16/01	ASYMPTOMATIC	IVDA	TRANSFUSION	45	12/12/01	7.0	12/12/01	31	12/12/01



# HBV LONGITUDINAL PANEL

(Catalog No. 6507)

For Investigational Research Use Only

**Application Notes:** This 14 member panel was collected from a 17 year old Caucasian Male.

PANEL NUMBER	CONCURRENT MEDICAL CONDITIONS	MEDICATIONS	ABBOTT HBSAG S/CO	ABBOTT HBC TOTAL S/CO	ABBOTT HBC IGM S/CO	DIASORI N HBEAG S/CO	DIASORI N ANTI-HBE S/CO	ABBOTT ANTI-HBS S/CO	ABBOTT ANTI-HIV 1/2 S/CO	ABBOTT HCV S/CO
6507-01	NONE	AGUA MINERALE/JUICE MULTIVITAMIN/ VITAMINAEVIT/VITAMINACIDIAS CORBINICI/SORBENT/SOL, RIBOXINI	>40.8	0.04	>6.9	0.7	0.6	0.3	0.3	3.8
6507-02	NONE	AGUA MINERALE/ JUICE MULTIVITAMIN	>40.8	0.06	>6.9	0.2	0.6	0.3	N/A	N/A
6507-03	NONE	DUOVIT	>40.8	0.07	6.0	0.1	0.6	0.2	N/A	N/A
6507-04	NONE	AGUA MINERALE/DUOVIT	>40.8	0.08	5.9	0.2	0.7	0.2	N/A	N/A
6507-05	NONE	NONE	0.04	0.08	4.0	0.3	0.3	0.4	N/A	N/A
6507-06	NONE	NONE	0.06	0.09	3.9	0.1	0.2	0.5	N/A	N/A
6507-07	NONE	HEPATIC TEA	0.04	0.09	3.7	0.2	0.2	0.5	N/A	N/A
6507-08	NONE	DIET N5	0.15	0.04	3.8	0.3	0.1	1.9	N/A	N/A
6507-09	NONE	DIET N5	0.15	0.03	2.3	0.3	0.3	0.55	N/A	N/A
6507-10	HCV	DIET N5	0.14	0.03	2.7	0.4	0.1	5.3	N/A	N/A
6507-11	HCV	DIET N5/ HEPATIC TEA	0.09	0.04	2.6	0.4	0.1	4.8	N/A	N/A
6507-12	HCV	DIET N5	0.02	0.04	2.3	0.4	0.1	4.4	N/A	N/A
6507-13	HCV	DIET N5	0.04	0.08	1.6	0.3	0.03	17.8	N/A	N/A
6507-14	HCV	DIET N5	0.09	0.09	1.4	0.3	0.04	19.0	N/A	N/A