



HBV LONGITUDINAL PANEL

(Catalog No. 6508)

For Investigational Research Use Only

Application Notes: This 15 member panel was collected from a 17 year old Caucasian Male.

PANEL NUMBER	DATE OF DRAW	DATE CRF COMPLETED	DATE OF BIRTH	DODX	SIGNS AND SYMPTOMS	RISK GROUP	MODE OF TRANSMISSION	ALT	ALT DATE	BILIRUBIN	BILIRUBIN DATE	SGOT	SGOT DATE
6508-01	5/16/01	5/18/01	1/12/84	5/7/01	JAUNDICE/ANOREXIA/HEPATOMEGALY/ DARK URINE/MALAISE	IVDA	TRANSFUSION	790	5/16/01	208	5/16/01	630	5/16/01
6508-02	5/28/01	5/29/01	1/12/84	5/7/01	JAUNDICE/FATIGUE/ HEPATOMEGALY/DARK URINE/MALAISE	IVDA	TRANSFUSION	630	5/28/01	160	5/28/01	540	5/28/01
6508-03	6/11/01	6/12/01	1/12/84	5/7/01	JAUNDICE/ HEPATOMEGALY/DARK URINE/ MALAISE	IVDA	TRANSFUSION	663	6/11/01	155	6/11/01	571	6/11/01
6508-04	6/28/01	7/3/01	1/12/84	5/7/01	ASYMPTOMATIC	IVDA	TRANSFUSION	43	6/29/01	34.2	6/29/01	44	6/29/01
6508-05	7/11/01	7/13/01	1/12/84	5/7/01	ASYMPTOMATIC	IVDA	TRANSFUSION	68.0	7/12/01	20.0	7/12/01	52	7/12/01
6508-06	7/25/01	7/27/01	1/12/84	5/7/01	ASYMPTOMATIC	IVDA	TRANSFUSION	16.3	7/26/01	10.2	7/26/01	26.5	7/26/01
6508-07	8/8/01	8/13/01	1/12/84	5/7/01	ASYMPTOMATIC	IVDA	TRANSFUSION	9.3	8/15/01	11.2	8/15/01	15.0	8/15/01
6508-08	8/22/01	8/26/01	1/12/84	5/7/01	ASYMPTOMATIC	IVDA	TRANSFUSION	21.3	8/23/01	13.1	8/23/01	16.7	8/23/01
6508-09	9/5/01	9/8/01	1/12/84	5/7/01	ASYMPTOMATIC	IVDA	TRANSFUSION	11.9	9/7/02	15.4	9/7/01	25.7	9/7/01
6508-10	9/19/01	9/24/01	1/12/84	5/7/01	ASYMPTOMATIC	IVDA	TRANSFUSION	16.7	9/21/01	13.4	9/21/01	127	9/21/01
6508-11	10/3/01	10/5/01	1/12/84	5/7/01	ASYMPTOMATIC	IVDA	TRANSFUSION	18.0	10/4/01	15.7	10/4/01	140.0	10/4/01
6508-12	10/16/01	10/19/01	1/12/84	5/7/01	ASYMPTOMATIC	IVDA	TRANSFUSION	10.1	10/17/01	12.4	10/17/01	26.0	10/17/01
6508-13	11/20/01	11/26/01	1/12/84	5/7/01	ASYMPTOMATIC	IVDA	TRANSFUSION	19	11/22/01	23.2	11/22/01	14	11/22/01
6508-14	12/11/01	12/12/01	1/12/84	5/7/01	ASYMPTOMATIC	IVDA	TRANSFUSION	8	12/12/01	17.8	12/12/01	18	12/12/01
6508-15	12/25/01	12/27/01	1/12/84	5/7/01	ASYMPTOMATIC	IVDA	TRANSFUSION	12	12/26/01	15.4	12/26/01	22	12/26/01



HBV LONGITUDINAL PANEL

(Catalog No. 6508)

For Investigational Research Use Only

Application Notes: This 15 member panel was collected from a 17 year old Caucasian Male.

PANEL NUMBER	MEDICATIONS	CONCURRENT MEDICAL CONDITIONS	ABBOTT HBSAG S/CO	ABBOTT HBC TOTAL S/CO	ABBOTT HBC IGM S/CO	DIASORIN HBEAG S/CO	DIASORIN ANTI-HBE S/CO	ABBOTT ANTI-HBS S/CO	ABBOTT ANTI-HIV 1/2 S/CO	ABBOTT HCV S/CO
6508-01	AGUA MINERALE/ SOL.ROBOXINI/AEVIT/SORBENT	NONE	>40.8	0.04	>6.9	0.3	0.6	0.3	0.2	0.3
6508-02	AGUA MINERALE/SORBENT/AEVIT/ ESSENCIALE	NONE	>40.8	0.05	>6.9	0.3	0.9	0.3	N/A	N/A
6508-03	AGUA MINERALE/HEPATIC TEA/DUOVIT	NONE	>40.8	0.04	>6.9	0.3	0.9	0.3	N/A	N/A
6508-04	NONE	NONE	>40.8	0.05	>6.9	0.3	1.0	0.2	N/A	N/A
6508-05	NONE	NONE	>40.8	0.06	>6.9	0.5	0.6	0.3	N/A	N/A
6508-06	DIET N5	NONE	7.1	0.04	3.1	0.07	0.1	0.4	N/A	N/A
6508-07	DIET N5	NONE	0.08	0.03	3.0	0.37	0.1	1.2	N/A	N/A
6508-08	DIET N5	NONE	0.04	0.03	3.4	0.07	0.1	1.95	N/A	N/A
6508-09	DIET N5	NONE	0.09	0.04	2.9	0.12	0.1	4.3	N/A	N/A
6508-10	DIET N5	NONE	0.09	0.03	2.9	0.08	0.1	7.5	N/A	N/A
6508-11	DIET N5	NONE	0.09	0.04	2.9	0.14	0.2	10.4	N/A	N/A
6508-12	DIET N5	NONE	0.04	0.07	2.5	0.04	0.2	9.3	N/A	N/A
6508-13	DIET N5	NONE	0.08	0.08	2.2	0.12	0.2	13.0	N/A	N/A
6508-14	DIET N5	NONE	0.02	0.08	2.1	0.08	0.3	10.9	N/A	N/A
6508-15	DIET N5	NONE	0.13	0.07	2.2	0.12	0.4	11.6	N/A	N/A