



HBV LONGITUDINAL PANEL

(Catalog No. 6509)

For Investigational Research Use Only

Application Notes: This 12 member panel was collected from a 18 year old Caucasian Male.

PANEL NUMBER	DATE OF DRAW	DATE CRF COMPLETED	DATE OF BIRTH	DODX	SIGNS AND SYMPTOMS	RISK GROUP	MODE OF TRANSMISSION	ALT	ALT DATE	BILIRUBIN	BILIRUBIN DATE	SGOT	SGOT DATE
6509-01	5/16/01	5/21/01	8/29/83	5/10/01	JAUNDICE/ ANOREXIA/FATIGUE/ HEPATOMEGALY/NAUSEA/ VOMTING/DARK URINE/ MALAISE/ABDOMINAL PAIN	IVDA	TRANSFUSION	4490	5/17/01	230	5/17/01	4680	5/17/01
6509-02	5/28/01	5/31/01	8/29/83	5/10/01	JAUNDICE/ ANOREXIA/ FATIGUE/HEPATOMEGALY/ DARK URINE/ MALAISE/PRURITIS	IVDA	TRANSFUSION	3600	5/28/01	190	5/28/01	3700	5/28/01
6509-03	6/13/01	6/14/01	8/29/83	5/10/01	JAUNDICE/ ANOREXIA/HEPATOMEGALY/ DARK URINE/MALAISE	IVDA	TRANSFUSION	3158	6/13/01	146	6/13/01	2338	6/13/01
6509-04	6/27/01	7/4/01	8/29/83	5/10/01	JAUNDICE/ HEPATOMEGALY/DARK URINE	IVDA	TRANSFUSION	820	6/28/01	100	6/28/01	600	6/28/01
6509-05	7/11/01	7/16/01	8/29/83	5/10/01	JAUNDICE/ HEPATOMEGALY/DARK URINE	IVDA	TRANSFUSION	430	7/12/01	76.0	7/12/01	290	7/12/01
6509-06	7/25/01	7/27/01	8/29/83	5/10/01	ASYMPTOMATIC	IVDA	TRANSFUSION	118.6	7/26/01	19.7	7/26/01	127.1	7/26/01
6509-07	8/8/01	8/13/01	8/29/83	5/10/01	ASYMPTOMATIC	IVDA	TRANSFUSION	13.5	8/15/01	15.3	8/15/01	19.5	8/15/01
6509-08	8/22/01	8/26/01	8/29/83	5/10/01	ASYMPTOMATIC	IVDA	TRANSFUSION	19.4	8/23/01	15.3	8/23/01	24.3	8/23/01
6509-09	9/5/01	9/8/01	8/29/83	5/10/01	ASYMPTOMATIC	IVDA	TRANSFUSION	14.4	9/7/01	13.3	9/7/01	24.7	9/7/01
6509-10	9/19/01	9/24/01	8/29/83	5/10/01	ASYMPTOMATIC	IVDA	TRANSFUSION	17.8	9/20/01	18.7	9/20/01	64.9	9/20/01
6509-11	10/3/01	10/5/01	8/29/83	5/10/01	ASYMPTOMATIC	IVDA	TRANSFUSION	15.4	10/4/01	16.4	10/4/01	40.0	10/4/01
6509-12	10/16/01	10/19/01	8/29/83	5/10/01	ASYMPTOMATIC	IVDA	TRANSFUSION	11.0	10/17/01	12.0	10/17/01	26.4	10/17/01



HBV LONGITUDINAL PANEL

(Catalog No. 6509)

For Investigational Research Use Only

Application Notes: This 12 member panel was collected from a 18 year old Caucasian Male.

PANEL NUMBER	CONCURRENT MEDICAL CONDITIONS	MEDICATIONS	ABBOTT HBSAG S/CO	ABBOTT HBC TOTAL S/CO	ABBOTT HBC IGM S/CO	DIASORIN HBEAG S/CO	DIASORIN ANTI-HBE S/CO	ABBOTT ANTI-HBS S/CO	ABBOTT ANTI-HIV 1/2 S/CO	ABBOTT HCV S/CO
6509-01	NONE	AGUA MINERALE/ ENTEROSORBENT/SOL. RIBOXINI/DUOVIT	>40.8	0.04	>6.9	1.1	0.9	0.3	0.2	0.3
6509-02	NONE	AGUA MINERALE/SORBENT/ SOL. PAPAVERINI/AEVIT/ESSENCIALE	>40.8	0.08	>6.9	0.4	1.0	0.3	N/A	N/A
6509-03	NONE	AGUAMINERALE/ENTEROSORENT/D UOVIT	>40.8	0.05	>6.9	0.4	1.3	0.3	N/A	N/A
6509-04	NONE	AGUA MINERALE/RIBOXINI/ ENTEROSORBENT	>40.8	0.05	>6.9	0.4	1.3	0.2	N/A	N/A
6509-05	NONE	AGUA MINERALE/HEPATIC TEA	>40.8	0.08	6.8	0.2	0.5	0.3	N/A	N/A
6509-06	NONE	DIET N5/ HEPATIC TEA	5.8	0.04	4.3	0.14	0.2	0.3	N/A	N/A
6509-07	NONE	DIET N5	0.04	0.04	2.7	0.09	0.1	0.93	N/A	N/A
6509-08	NONE	DIET N5	0.09	0.04	2.8	0.03	0.1	1.9	N/A	N/A
6509-09	NONE	DIET N5	0.15	0.56	2.1	0.08	0.2	2.4	N/A	N/A
6509-10	NONE	DIET N5	0.49	0.05	2.1	0.05	0.2	6.9	N/A	N/A
6509-11	NONE	DIET N5	0.05	0.05	1.9	0.08	0.2	8.1	N/A	N/A
6509-12	NONE	DIET N5	0.02	0.07	2.0	0.05	0.2	8.1	N/A	N/A