



HBV LONGITUDINAL PANEL

(Catalog No. 6510)

For Investigational Research Use Only

Application Notes: This 12 member panel was collected from a 32 year old Caucasian Male.

PANEL NUMBER	DATE OF DRAW	DATE CRF COMPLETED	DATE OF BIRTH	DODX	SIGNS AND SYMPTOMS	RISK GROUP	MODE OF TRANSMISSION	ALT	ALT DATE	BILIRUBIN	BILIRUBIN DATE	SGOT	SGOT DATE
6510-01	5/22/01	5/25/01	5/1/69	5/22/01	JAUNDICE/ANOREXIA/FATIGUE/HEPATOMEGALY/DARK URINE/MALAISE/SKIN IRRITATION/PRURITIS	NO KNOWN RISK	NOT DETERMINED	1520	5/22/01	220.4	5/22/01	860	5/22/01
6510-02	6/5/01	6/7/01	5/1/69	5/22/01	JAUNDICE/ANOREXIA/FATIGUE/HEPATOMEGALY/NAUSEA/VOMITING/DARK URINE/MALAISE/PRURITIS	NO KNOWN RISK	NOT DETERMINED	1089	6/5/01	497	6/5/01	912	6/5/01
6510-03	6/19/01	6/22/01	5/1/69	5/22/01	JAUNDICE/ANOREXIA/HEPATOMEGALY/DARK URINE	NO KNOWN RISK	NOT DETERMINED	137	6/20/01	183.9	6/20/01	228	6/20/01
6510-04	7/3/01	7/6/01	5/1/69	5/22/01	JAUNDICE/ANOREXIA/HEPATOMEGALY/DARK URINE/MALAISE	NO KNOWN RISK	NOT DETERMINED	251	7/4/01	126	7/4/01	215	7/4/01
6510-05	7/17/01	7/20/01	5/1/69	5/22/01	JAUNDICE/HEPATOMEGALY/DARK URINE	NO KNOWN RISK	NOT DETERMINED	23.0	7/18/01	113.0	7/18/01	290	7/18/01
6510-06	7/31/01	8/6/01	5/1/69	5/22/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	32.4	8/3/01	27.5	8/3/01	42.3	8/3/01
6510-07	8/14/01	8/20/01	5/1/69	5/22/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	14.2	8/15/01	14.7	8/15/01	20.7	8/15/01
6510-08	8/28/01	9/1/01	5/1/69	5/22/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	9.5	8/28/01	18.7	8/28/01	19.6	8/28/01
6510-09	9/11/01	9/17/01	5/1/69	5/22/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	14.2	9/14/01	13.2	9/14/01	23.6	9/14/01
6510-10	9/25/01	9/26/01	5/1/69	5/22/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	17.3	9/26/01	12.9	9/26/01	35.3	9/26/01
6510-11	10/9/01	10/12/01	5/1/69	5/22/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	16.4	10/10/01	26.4	10/10/01	20.0	10/10/01
6510-12	10/22/01	10/24/01	5/1/69	5/22/01	ASYMPTOMATIC	NO KNOWN RISK	NOT DETERMINED	10.1	10/23/01	12.4	10/23/01	26.0	10/23/01



HBV LONGITUDINAL PANEL

(Catalog No. 6510)

For Investigational Research Use Only

Application Notes: This 12 member panel was collected from a 32 year old Caucasian Male.

PANEL NUMBER	CONCURRENT MEDICAL CONDITIONS	MEDICATIONS	ABBOTT HBSAG S/CO	ABBOTT HBC TOTAL S/CO	ABBOTT HBC IGM S/CO	DIASORIN HBEAG S/CO	DIASORIN ANTI-HBE S/CO	ABBOTT ANTI-HBS S/CO	ABBOTT ANTI-HIV 1/2 S/CO	ABBOTT HCV S/CO
6510-01	NONE	AGUA MINERALE/AEVIT/ SOL. PAPAVERINI/SOL. RIBOXINI/SORBENT/ SL. SUPRASTINI	34.7	0.03	>6.9	21.8	2.5	0.2	0.2	0.2
6510-02	NONE	AGUA MINERALE/ ENTEROSORBENT/SOL. RIBOXINI/PLASMAFERES/ DUOVIT/SOL. SUPRASTINI/MESIN	>40.8	0.03	>6.9	21.0	2.7	0.3	N/A	N/A
6510-03	NONE	AGUA MINERALE/ ENTEROSORBENT/PLASMAFERES/NO-SPANI	>40.8	0.04	4.8	6.8	1.7	0.3	N/A	N/A
6510-04	NONE	AGUA MINERALE/ ENTEROSORBENT/PLASMAFERESI/DUOVIT	>40.8	0.04	4.1	5.6	1.8	0.3	N/A	N/A
6510-05	NONE	DIET N5/ AQUA MINERALE/ HEPATIC TEA/ AEVIT/ ENTEROSOBENT	>38	0.03	2.7	17.2	1.6	0.15	N/A	N/A
6510-06	NONE	DIET N5	>38	0.04	5.5	0.02	0.2	0.4	N/A	N/A
6510-07	NONE	DIET N5	0.08	0.04	2.9	0.08	0.1	0.84	N/A	N/A
6510-08	NONE	DIET N5	0.23	0.04	2.8	0.02	0.1	1.54	N/A	N/A
6510-09	NONE	DIET N5	0.09	0.04	1.205	0.16	0.2	3.4	N/A	N/A
6510-10	NONE	DIET N5	0.11	0.04	1.8	0.06	0.2	6.2	N/A	N/A
6510-11	NONE	DIET N5	0.08	0.06	2.0	0.05	0.3	5.7	N/A	N/A
6510-12	NONE	DIET N5	0.04	0.06	1.96	0.02	0.2	6.0	N/A	N/A