



HBV LONGITUDINAL PANEL

(Catalog No. 6512)

For Investigational Research Use Only

Application Notes: This 15 member panel was collected from a 18 year old Caucasian Male.

PANEL NUMBER	DATE OF DRAW	DATE CRF COMPLETED	DATE OF BIRTH	DODX	SIGNS AND SYMPTOMS	RISK GROUP	MODE OF TRANSMISSION	ALT	ALT DATE	BILIRUBIN	BILIRUBIN DATE	SGOT	SGOT DATE
6512-01	5/30/01	5/31/01	11/16/83	5/29/01	JAUNDICE/ HEPATOMEGALY/ DARK URINE	IVDA	TRANSFUSION	180	5/30/01	42.5	5/30/01	475	5/30/01
6512-02	6/13/01	6/14/01	11/16/83	5/29/01	ASYMPTOMATIC	IVDA	TRANSFUSION	234	6/13/01	19.7	6/13/01	130	6/13/01
6512-03	6/27/01	7/3/01	11/16/83	5/29/01	ASYMPTOMATIC	IVDA	TRANSFUSION	29.8	6/28/01	18.7	6/28/01	34.6	6/28/01
6512-04	7/11/01	7/13/01	11/16/83	5/29/01	ASYMPTOMATIC	IVDA	TRANSFUSION	35.4	7/12/01	23.6	7/12/01	44.2	7/12/01
6512-05	7/25/01	7/27/01	11/16/83	5/29/01	ASYMPTOMATIC	IVDA	TRANSFUSION	18.2	7/26/01	17.7	7/26/01	35.6	7/26/01
6512-06	8/8/01	8/13/01	11/16/83	5/29/01	ASYMPTOMATIC	IVDA	TRANSFUSION	24.1	8/15/01	18.7	8/15/01	12.9	8/15/01
6512-07	8/22/01	8/26/01	11/16/83	5/29/01	ASYMPTOMATIC	IVDA	TRANSFUSION	13.0	8/23/01	19.7	8/23/01	27.6	8/23/01
6512-08	9/5/01	9/7/01	11/16/83	5/29/01	ASYMPTOMATIC	IVDA	TRANSFUSION	20.0	9/7/01	21.5	9/7/01	25.0	9/7/01
6512-09	9/19/01	9/21/01	11/16/83	5/29/01	ASYMPTOMATIC	IVDA	TRANSFUSION	12.3	9/20/01	17.2	9/20/01	24.2	9/20/01
6512-10	10/3/01	10/5/01	11/16/83	5/29/01	ASYMPTOMATIC	IVDA	TRANSFUSION	16.6	10/4/01	15.5	10/4/01	28.3	10/4/01
6512-11	10/17/01	10/19/01	11/16/83	5/29/01	ASYMPTOMATIC	IVDA	TRANSFUSION	9.0	10/18/01	21.8	10/18/01	29.9	10/18/01
6512-12	10/30/01	11/1/01	11/16/83	5/29/01	ASYMPTOMATIC	IVDA	TRANSFUSION	12.1	10/31/01	12.7	10/31/01	27.6	10/31/01
6512-13	12/11/01	12/12/01	11/16/83	5/29/01	ASYMPTOMATIC	IVDA	TRANSFUSION	7	12/12/01	21.4	12/12/01	20	12/12/01
6512-14	12/25/01	12/27/01	11/16/83	5/29/01	ASYMPTOMATIC	IVDA	TRANSFUSION	6	12/26/01	28.8	12/26/01	23	12/26/01
6512-15	1/8/02	1/11/02	11/16/83	5/29/01	ASYMPTOMATIC	IVDA	TRANSFUSION	15.0	1/9/02	18.7	1/9/02	17.0	1/9/02



HBV LONGITUDINAL PANEL

(Catalog No. 6512)

For Investigational Research Use Only

Application Notes: This 15 member panel was collected from a 18 year old Caucasian Male.

PANEL NUMBER	CONCURRENT MEDICAL CONDITIONS	MEDICATIONS	ABBOTT HBSAG S/CO	ABBOTT HBC TOTAL S/CO	ABBOTT HBC IGM S/CO	DIASORIN HBEAG S/CO	DIASORIN ANTI-HBE S/CO	ABBOTT ANTI-HBS S/CO	ABBOTT ANTI-HIV 1/2 S/CO	ABBOTT HCV S/CO
6512-01	NONE	AGUA MINERALE/DUOVIT/ SOL. RIBOXINI	>40.8	0.07	>6.9	0.5	0.6	0.3	0.2	2.5
6512-02	NONE	DUOVIT	>40.8	0.08	>6.9	0.4	0.9	0.5	N/A	N/A
6512-03	NONE	NONE	18.8	0.09	>6.9	0.2	1.0	0.3	N/A	N/A
6512-04	NONE	NONE	1.6	0.11	>6.9	0.2	1.0	0.3	N/A	N/A
6512-05	NONE	DIET N5	0.3	0.03	>6.5	0.3	0.3	0.4	N/A	N/A
6512-06	NONE	DIET N5	0.1	0.03	5.9	0.2	0.4	0.5	N/A	N/A
6512-07	NONE	DIET N5	0.1	0.04	5.9	0.3	0.3	0.6	N/A	N/A
6512-08	NONE	DIET N5	0.1	0.04	5.1	0.2	0.2	0.7	N/A	N/A
6512-09	NONE	DIET N5	0.8	0.03	4.7	0.3	0.2	0.9	N/A	N/A
6512-10	NONE	DIET N5	0.02	0.03	4.1	0.2	0.2	1.1	N/A	N/A
6512-11	NONE	DIET N5	0.1	0.05	5.8	0.3	0.2	1.1	N/A	N/A
6512-12	NONE	DIET N5	0.02	0.04	5.4	0.2	0.2	1.2	N/A	N/A
6512-13	NONE	DIET N5	0.1	0.06	4.4	0.4	0.2	1.6	N/A	N/A
6512-14	NONE	DIET N5	0.1	0.06	4.1	0.2	0.3	1.5	N/A	N/A
6512-15	NONE	DIET N5	0.1	0.05	4.1	0.3	0.4	1.31	N/A	N/A