



HBV LONGITUDINAL PANEL

(Catalog No. 6536)

For Investigational Research Use Only

Application Notes: This 12 member panel was collected from a 24 year old Caucasian Male.

PANEL NUMBER	DATE OF DRAW	DATE CRF COMPLETED	DATE OF BIRTH	DODX	SIGNS AND SYMPTOMS	RISK GROUP	MODE OF TRANSMISSION	ALT	ALT DATE	BILIRUBIN	BILIRUBIN DATE	SGOT	SGOT DATE
6536-01	11/13/01	11/16/01	5/27/77	11/9/01	JAUNDICE/ HEPATOMEGALY	IVDA	TRANSFUSION	574.0	11/15/01	46.5	11/15/01	186.0	11/15/01
6536-02	11/26/01	11/28/01	5/27/77	11/9/01	HEPATOMEGALY	IVDA	TRANSFUSION	220.9	11/27/01	25.9	11/27/01	102.6	11/27/01
6536-03	12/11/01	12/14/01	5/27/77	11/9/01	ASYMPTOMATIC	IVDA	TRANSFUSION	195.0	12/12/01	14.3	12/12/01	21.0	12/12/01
6536-04	12/25/01	12/27/01	5/27/77	11/9/01	ASYMPTOMATIC	IVDA	TRANSFUSION	45.0	12/26/01	14.3	12/26/01	33.0	12/26/01
6536-05	1/8/02	1/10/02	5/27/77	11/9/01	ASYMPTOMATIC	IVDA	TRANSFUSION	38.0	1/9/02	11.0	1/9/02	31.0	1/9/02
6536-06	1/22/02	1/24/02	5/27/77	11/9/01	ASYMPTOMATIC	IVDA	TRANSFUSION	40.0	1/23/02	6.8	1/23/02	29.0	1/23/02
6536-07	2/5/02	2/7/02	5/27/77	11/9/01	ASYMPTOMATIC	IVDA	TRANSFUSION	24.0	2/6/02	10.7	2/6/02	31.6	2/6/02
6536-08	2/19/02	2/22/02	5/27/77	11/9/01	ASYMPTOMATIC	IVDA	TRANSFUSION	118.3	2/20/02	7.5	2/20/02	68.8	2/20/02
6536-09	3/5/02	3/7/02	5/27/77	11/9/01	ASYMPTOMATIC	IVDA	TRANSFUSION	5.8	3/6/02	10.0	3/6/02	11.8	3/6/02
6536-10	3/19/02	3/21/02	5/27/77	11/9/01	ASYMPTOMATIC	IVDA	TRANSFUSION	68.0	3/20/02	3.7	3/20/02	43.0	3/20/02
6536-11	4/2/02	4/5/02	5/27/77	11/9/01	ASYMPTOMATIC	IVDA	TRANSFUSION	25.0	4/4/02	6.6	4/4/02	35.0	4/4/02
6536-12	4/16/02	4/20/02	5/27/77	11/9/01	ASYMPTOMATIC	IVDA	TRANSFUSION	67.0	4/17/02	4.3	4/17/02	41.0	4/17/02



HBV LONGITUDINAL PANEL

(Catalog No. 6536)

For Investigational Research Use Only

Application Notes: This 12 member panel was collected from a 24 year old Caucasian Male.

PANEL NUMBER	CONCURRENT MEDICAL CONDITIONS	MEDICATIONS	ABBOTT HBSAG S/CO	ABBOTT HBC TOTAL S/CO	ABBOTT HBC IGM S/CO	HBEAG DIASORIN S/CO	DIASORIN ANTI-HBE S/CO	ABBOTT ANTI-HBS S/CO	ABBOTT ANTI-HIV 1/2 S/CO	ABBOTT HCV S/CO
6536-01	HCV	DIET N5/ AQUA MINERALE/ AEVIT/ HEPATIC TEA	>36.4	0.04	>4.7	0.31	0.3	0.3	0.09	>6.2
6536-02	HCV	DIET N5/ HEPATIC TEA/ AEVIT	>36.4	0.03	>4.7	0.31	0.4	0.3	N/A	N/A
6536-03	HCV	DIET N5/ HEPATIC TEA/ ESSENCIALE	>36.4	0.07	>4.7	0.34	0.1	0.4	N/A	N/A
6536-04	HCV	DIET N5/ ESSENCIALE	5.6	0.05	>4.7	0.26	0.07	0.3	N/A	N/A
6536-05	HCV	DIET N5	0.4	0.11	3.9	0.07	0.05	0.5	N/A	N/A
6536-06	HCV	DIET N5	0.7	0.09	3.1	0.16	0.05	0.7	N/A	N/A
6536-07	HCV	DIET N5	0.2	0.06	2.5	0.09	0.03	0.9	N/A	N/A
6536-08	HCV	DIET N5	0.2	0.07	2.4	0.23	0.04	1.3	N/A	N/A
6536-09	HCV	DIET N5	0.2	0.07	2.1	0.16	0.04	1.8	N/A	N/A
6536-10	HCV	DIET N5	0.2	0.07	1.8	0.26	0.05	2.2	N/A	N/A
6536-11	HCV	DIET N5	0.1	0.07	2.0	0.16	0.05	2.2	N/A	N/A
6536-12	HCV	DIET N5	0.2	0.06	1.7	0.24	0.05	2.1	N/A	N/A