



HBV LONGITUDINAL PANEL

(Catalog No. 6541)

For Investigational Research Use Only

Application Notes: This 12 member panel was collected from a 20 year old Caucasian Male.

PANEL NUMBER	DATE OF DRAW	DATE CRF COMPLETED	DATE OF BIRTH	DODX	SIGNS AND SYMPTOMS	RISK GROUP	MODE OF TRANSMISSION	ALT	ALT DATE	BILIRUBIN	BILIRUBIN DATE	SGOT	SGOT DATE
6541-01	12/3/01	12/5/01	8/13/81	12/1/01	HEPATOMEGALY	IVDA	TRANSFUSION	955.0	12/4/01	18.3	12/4/01	473.9	12/4/01
6541-02	12/17/01	12/20/01	8/13/81	12/1/01	HEPATOMEGALY	IVDA	TRANSFUSION	466.0	12/18/01	9.8	12/18/01	138.0	12/18/01
6541-03	12/31/01	1/10/02	8/13/81	12/1/01	ASYMPTOMATIC	IVDA	TRANSFUSION	28.0	1/3/02	15.6	1/3/02	78.0	1/3/02
6541-04	1/14/02	1/18/02	8/13/81	12/1/01	ASYMPTOMATIC	IVDA	TRANSFUSION	24.0	1/15/02	6.6	1/15/02	21.0	1/15/02
6541-05	1/28/02	1/30/02	8/13/81	12/1/01	ASYMPTOMATIC	IVDA	TRANSFUSION	12.0	1/29/02	10.2	1/29/02	18.0	1/29/02
6541-06	2/11/02	2/13/02	8/13/81	12/1/01	ASYMPTOMATIC	IVDA	TRANSFUSION	18.0	2/12/02	9.0	2/12/02	46.2	2/12/02
6541-07	2/25/02	2/27/02	8/13/81	12/1/01	ASYMPTOMATIC	IVDA	TRANSFUSION	3.0	2/26/02	10.2	2/26/02	22.0	2/26/02
6541-08	3/11/02	3/12/02	8/13/81	12/1/01	ASYMPTOMATIC	IVDA	TRANSFUSION	41.4	3/12/02	16.9	3/12/02	53.1	3/12/02
6541-09	3/25/02	3/27/02	8/13/81	12/1/01	ASYMPTOMATIC	IVDA	TRANSFUSION	11.0	3/26/02	4.4	3/26/02	16.0	3/26/02
6541-10	4/8/02	4/10/02	8/13/81	12/1/01	ASYMPTOMATIC	IVDA	TRANSFUSION	19.0	4/9/02	9.4	4/9/02	10.0	4/9/02
6541-11	4/22/02	4/24/02	8/13/81	12/1/01	ASYMPTOMATIC	IVDA	TRANSFUSION	13.0	4/23/02	6.8	4/23/02	7.0	4/23/02
6541-12	5/6/02	5/8/02	8/13/81	12/1/01	ASYMPTOMATIC	IVDA	TRANSFUSION	14.6	5/7/02	11.0	5/7/02	25.3	5/7/02



HBV LONGITUDINAL PANEL

(Catalog No. 6541)

For Investigational Research Use Only

Application Notes: This 12 member panel was collected from a 20 year old Caucasian Male.

PANEL NUMBER	CONCURRENT MEDICAL CONDITIONS	MEDICATIONS	ABBOTT HBSAG S/CO	ABBOTT HBC TOTAL S/CO	ABBOTT HBC IGM S/CO	DIASORIN HBEAG S/CO	DIASORIN ANTI-HBE S/CO	ABBOTT ANTI-HBS S/CO	ABBOTT ANTI-HIV 1/2 S/CO	ABBOTT HCV S/CO
6541-01	NONE	DIET N5/ HEPATIC TEA/ AEVIT/ RIBOXINI	>37	0.1	>6.5	1.8	0.5	0.3	0.09	0.2
6541-02	NONE	DIET N5/ HEPATIC TEA/ AEVIT/ RIBOXINI	>37	0.1	>6.5	0.39	0.3	0.5	N/A	N/A
6541-03	NONE	DIET N5	>37	0.1	>6.5	0.12	0.3	0.5	N/A	N/A
6541-04	NONE	DIET N5	1.1	0.1	6.0	0.12	0.1	0.5	N/A	N/A
6541-05	NONE	DIET N5	0.1	0.1	3.9	0.16	0.1	0.5	N/A	N/A
6541-06	NONE	DIET N5	0.1	0.1	3.3	0.11	0.09	0.8	N/A	N/A
6541-07	NONE	DIET N5	0.2	0.3	3.1	0.34	0.09	1.08	N/A	N/A
6541-08	NONE	DIET N5	0.1	0.1	2.6	0.18	0.08	2.5	N/A	N/A
6541-09	NONE	DIET N5	0.1	0.1	1.9	0.04	0.06	2.3	N/A	N/A
6541-10	NONE	DIET N5	0.1	0.1	1.8	0.16	0.07	2.5	N/A	N/A
6541-11	NONE	DIET N5	0.1	0.1	1.4	0.12	0.05	3.3	N/A	N/A
6541-12	NONE	DIET N5	0.1	0.1	1.6	0.09	0.05	2.2	N/A	N/A