



HBV LONGITUDINAL PANEL

(Catalog No. 6542)

For Investigational Research Use Only

Application Notes: This 12 member panel was collected from a 17 year old Caucasian Male.

PANEL NUMBER	DATE OF DRAW	DATE CRF COMPLETED	DATE OF BIRTH	DODX	SIGNS AND SYMPTOMS	RISK GROUP	MODE OF TRANSMISSION	ALT	ALT DATE	BILIRUBIN	BILIRUBIN DATE	SGOT	SGOT DATE	CONCURRENT MEDICAL CONDITIONS
6542-01	12/25/01	12/27/01	10/8/84	12/25/01	HEPATOMEGALY	IVDA	TRANSFUSION	401	12/26/01	18.2	12/26/01	216	12/26/01	NONE
6542-02	1/8/02	1/10/02	10/8/84	12/25/01	HEPATOMEGALY	IVDA	TRANSFUSION	52.0	1/9/02	18.3	1/9/02	36.0	1/9/02	NONE
6542-03	1/22/02	1/24/02	10/8/84	12/25/01	ASYMPTOMATIC	IVDA	TRANSFUSION	24.5	1/23/02	4.3	1/23/02	28.4	1/23/02	NONE
6542-04	2/5/02	2/7/02	10/8/84	12/25/01	ASYMPTOMATIC	IVDA	TRANSFUSION	12.8	2/6/02	10.2	2/6/02	24.5	2/6/02	NONE
6542-05	2/19/02	2/22/02	10/8/84	12/25/01	ASYMPTOMATIC	IVDA	TRANSFUSION	13.3	2/20/02	12.3	2/20/02	25.7	2/20/02	NONE
6542-06	3/5/02	3/6/02	10/8/84	12/25/01	ASYMPTOMATIC	IVDA	TRANSFUSION	8.8	3/6/02	13.3	3/6/02	43.5	3/6/02	NONE
6542-07	3/19/02	3/21/02	10/8/84	12/25/01	ASYMPTOMATIC	IVDA	TRANSFUSION	20.0	3/20/02	6.7	3/20/02	13.0	3/20/02	NONE
6542-08	4/2/02	4/5/02	10/8/84	12/25/01	ASYMPTOMATIC	IVDA	TRANSFUSION	32.0	4/4/02	5.5	4/4/02	38.0	4/4/02	NONE
6542-09	4/16/02	4/19/02	10/8/84	12/25/01	ASYMPTOMATIC	IVDA	TRANSFUSION	253.0	4/17/02	5.9	4/17/02	168.0	4/17/02	NONE
6542-10	4/30/02	5/4/02	10/8/84	12/25/01	HEPATOMEGALY	IVDA	TRANSFUSION	221.0	5/3/02	13.4	5/3/02	116.0	5/3/02	NONE
6542-11	5/14/02	5/16/02	10/8/84	12/25/01	HEPATOMEGALY	IVDA	TRANSFUSION	220.0	5/15/02	11.8	5/15/02	118.0	5/15/02	NONE
6542-12	5/28/02	5/30/02	10/8/84	12/25/01	HEPATOMEGALY	IVDA	TRANSFUSION	226.0	5/29/02	9.8	5/29/02	165.4	5/29/02	NONE



HBV LONGITUDINAL PANEL

(Catalog No. 6542)

For Investigational Research Use Only

Application Notes: This 12 member panel was collected from a 17 year old Caucasian Male.

PANEL NUMBER	CONCURRENT MEDICAL CONDITIONS	MEDICATIONS	ABBOTT HBSAG S/CO	ABBOTT HBC TOTAL S/CO	ABBOTT HBC IGM S/CO	DIASORIN HBEAG S/CO	DIASORIN ANTI-HBE S/CO	ABBOTT ANTI-HBS S/CO	ABBOTT ANTI-HIV 1/2 S/CO	ABBOTT HCV S/CO
6542-01	NONE	DIET N5/ HEPATIC TEA/ ESSENCIALE	>37	0.04	>6.5	0.31	0.1	0.4	0.13	5.5
6542-02	NONE	DIET N5/ AQUA MINERALE/ AEVIT	>37	0.1	>6.5	0.19	0.1	0.2	N/A	N/A
6542-03	NONE	DIET N5	>37	0.1	4.9	0.28	0.2	0.3	N/A	N/A
6542-04	NONE	DIET N5	2.9	0.1	3.4	0.19	0.1	0.3	N/A	N/A
6542-05	NONE	DIET N5	0.275	0.1	2.4	0.07	0.07	0.4	N/A	N/A
6542-06	NONE	DIET N5	0.125	0.1	1.8	0.11	0.06	0.4	N/A	N/A
6542-07	NONE	DIET N5	2.05	0.1	1.4	0.05	0.09	0.6	N/A	N/A
6542-08	NONE	DIET N5	0.1	0.04	1.2	0.11	0.09	2.1	N/A	N/A
6542-09	NONE	DIET N5/ DUOVIT	0.1	0.1	1.1	0.09	0.1	2.4	N/A	N/A
6542-10	NONE	DIET N5/ VIFERONI/ HEPATIC TEA	0.1	0.04	0.945	0.15	0.2	3.5	N/A	N/A
6542-11	NONE	DIET N5/ VIFERONI/ HEPATIC TEA	0.1	0.1	0.96	0.01	0.2	4.8	N/A	N/A
6542-12	NONE	DIET N5	0.1	0.1	1.165	0.04	0.2	6.1	N/A	N/A