



HBV LONGITUDINAL PANEL

(Catalog No. 6545)

For Investigational Research Use Only

Application Notes: This 12 member panel was collected from a 18 year old Caucasian Male.

PANEL NUMBER	DATE OF DRAW	DATE CRF COMPLETED	DATE OF BIRTH	DODX	SIGNS AND SYMPTOMS	RISK GROUP	MODE OF TRANSMISSION	ALT	ALT DATE	BILIRUBIN	BILIRUBIN DATE	SGOT	SGOT DATE
6545-01	1/23/02	1/25/02	7/5/83	1/21/02	JAUNDICE/ HEPATOMEGALY	IVDA	TRANSFUSION	585.0	1/24/02	28.9	1/24/02	327.0	1/24/02
6545-02	2/6/02	2/8/02	7/5/83	1/21/02	ASYMPTOMATIC	IVDA	TRANSFUSION	278.0	2/7/02	19.8	2/7/02	128.0	2/7/02
6545-03	2/19/02	2/22/02	7/5/83	1/21/02	ASYMPTOMATIC	IVDA	TRANSFUSION	157	2/20/02	19.2	2/20/02	94.4	2/20/02
6545-04	3/6/02	3/7/02	7/5/83	1/21/02	ASYMPTOMATIC	IVDA	TRANSFUSION	10.3	3/6/02	26.5	3/6/02	42.9	3/6/02
6545-05	3/19/02	3/21/02	7/5/83	1/21/02	ASYMPTOMATIC	IVDA	TRANSFUSION	120.0	3/20/02	9.3	3/20/02	135.0	3/20/02
6545-06	4/2/02	4/4/02	7/5/83	1/21/02	ASYMPTOMATIC	IVDA	TRANSFUSION	76.0	4/3/02	8.3	4/3/02	36.0	4/3/02
6545-07	4/16/02	4/18/02	7/5/83	1/21/02	ASYMPTOMATIC	IVDA	TRANSFUSION	66.0	4/17/02	11.8	4/17/02	48.0	4/17/02
6545-08	4/30/02	5/4/02	7/5/83	1/21/02	ASYMPTOMATIC	IVDA	TRANSFUSION	104.7	5/3/02	10.0	5/3/02	72.8	5/3/02
6545-09	5/14/02	5/16/02	7/5/83	1/21/02	ASYMPTOMATIC	IVDA	TRANSFUSION	52.3	5/15/02	14.7	5/15/02	50.4	5/15/02
6545-10	5/28/02	5/30/02	7/5/83	1/21/02	ASYMPTOMATIC	IVDA	TRANSFUSION	51.5	5/29/02	15.3	5/29/02	48.5	5/29/02
6545-11	6/11/02	6/14/02	7/5/83	1/21/02	ASYMPTOMATIC	IVDA	TRANSFUSION	54.0	6/13/02	33.3	6/13/02	85.0	6/13/02
6545-12	6/25/02	6/27/02	7/5/83	1/21/02	ASYMPTOMATIC	IVDA	TRANSFUSION	35.7	6/26/02	17.3	6/26/02	34.3	6/26/02



HBV LONGITUDINAL PANEL

(Catalog No. 6545)

For Investigational Research Use Only

Application Notes: This 12 member panel was collected from a 18 year old Caucasian Male.

PANEL NUMBER	CONCURRENT MEDICAL CONDITIONS	MEDICATIONS	ABBOTT HBSAG S/CO	ABBOTT HBC TOTAL S/CO	ABBOTT HBC IGM S/CO	DIASORIN HBEAG S/CO	DIASORIN ANTI-HBE S/CO	ABBOTT ANTI-HBS S/CO	ABBOTT ANTI-HIV 1/2 S/CO	ABBOTT HCV S/CO
6545-01	NONE	DIET N5/ AQUA MINERALE/ RIBOXINI/ AEVIT	>36.4	0.04	>5.5	0.11	0.2	0.3	0.12	5.1
6545-02	NONE	DIET N5/ RIBOXINI/ AEVIT/ HEPATIC TEA	>36.4	0.03	>5.5	0.16	0.2	0.3	N/A	N/A
6545-03	NONE	DIET N5/ HEPATIC TEA	>36.4	0.1	>5.5	0.05	0.2	0.3	N/A	N/A
6545-04	NONE	DIET N5	>36.4	0.5	>5.5	0.18	0.07	0.5	N/A	N/A
6545-05	NONE	DIET N5/ HEPATIC TEA/ DUOVIT	>36.4	0.1	4.7	0.04	0.04	0.4	N/A	N/A
6545-06	NONE	DIET N5	>36.4	0.04	4.1	0.12	0.04	0.4	N/A	N/A
6545-07	NONE	DIET N5	>36.4	0.1	4.0	0.04	0.04	0.4	N/A	N/A
6545-08	NONE	DIET N5	>36.4	0.1	3.7	0.14	0.05	0.3	N/A	N/A
6545-09	NONE	DIET N5	24.2	0.1	3.4	0.07	0.04	0.5	N/A	N/A
6545-10	NONE	DIET N5	5.6	0.1	3.1	0.16	0.09	0.4	N/A	N/A
6545-11	NONE	DIET N5	1.7	0.1	3.2	0.08	0.04	0.3	N/A	N/A
6545-12	NONE	DIET N5	0.5	0.1	2.9	0.11	0.05	0.3	N/A	N/A