



## HAV VACCINE PANEL

(Catalog No. 9611-027-2794)

For Investigational Research Use Only

**Application Notes:** This 5 member panel was collected from a *40 year old Caucasian Female*.

LOT NUMBER	DATE OF VACCINE	DATE OF DRAW	HAV RESULT	HAV Total Abbott S/Co Pre	HAV Total Abbott S/Co Post
9611-027-2794 A	10/3/96	10/3/96	<b>2.191</b>	<b>2.191</b>	0.912
9611-027-2794 B	-	10/24/96	0.912	-	-
9611-027-2794 POST 90 DAYS	-	1/23/97	0.890	-	-
9611-027-2794 120 DAYS	-	2/13/97	0.087	-	-
9611-027-2794 150 DAYS	-	5/8/97	0.120	-	-