



HAV VACCINE PANEL

(Catalog No. 9611-027-2801

For Investigational Research Use Only

Application Notes: This 5 member panel was collected from a *25 year old Black Male*.

LOT NUMBER	DATE OF VACCINE	DATE OF DRAW	HAV RESULT	HAV Total Abbott S/Co Pre	HAV Total Abbott S/Co Post
9611-027-2801 A	10/3/96	10/3/96	2.188	2.118	1.567
9611-027-2801 B	-	10/24/96	1.567	-	-
9611-027-2801 POST 90 DAYS	-	1/23/97	1.948	-	-
9611-027-2801 120 DAYS	-	4/15/97	0.209	-	-
9611-027-2801 150 DAYS	-	5/8/97	0.108	-	-