



HAV VACCINE PANEL

(Catalog No. 9611-027-2803)

For Investigational Research Use Only

Application Notes: This 5 member panel was collected from a *27 year old Black Male*.

LOT NUMBER	DATE OF VACCINE	DATE OF DRAW	HAV RESULT	HAV Total Abbott S/Co Pre	HAV Total Abbott S/Co Post
9611-027-2803 A	10/3/96	10/3/96	2.173	2.173	1.781
9611-027-2803 B	-	10/24/96	1.781	-	-
9611-027-2803 POST 90 DAYS	-	1/23/97	1.679	-	-
9611-027-2803 120 DAYS	-	4/13/97	1.944	-	-
9611-027-2803 150 DAYS	-	5/8/97	0.123	-	-