



HAV VACCINE PANEL

(Catalog No. 9611-027-2812)

For Investigational Research Use Only

Application Notes: This 5 member panel was collected from a *18 year old Black Male*.

LOT NUMBER	DATE OF VACCINE	DATE OF DRAW	HAV RESULT	HAV Total Abbott S/Co Pre	HAV Total Abbott S/Co Post
9611-027-2812 A	10/3/96	10/3/96	2.241	2.241	1.125
9611-027-2812 B	-	10/24/96	1.125	-	-
9611-027-2812 POST 90 DAYS	-	4/2/97	0.128	-	-
9611-027-2812 120 DAYS	-	4/16/97	0.128	-	-
9611-027-2812 150 DAYS	-	5/13/97	0.143	-	-