



MIXED TITER SYPHILIS PANEL

(Catalog No. K-ZMC002)

For Investigational Research Use Only

Application Notes: The syphilis positive donor plasmas used for this panel were obtained under informed consent and collected in US FDA licensed plasma centers. Donors were untreated and asymptomatic at the time of collection. Positive test results below are highlighted in red.

Member Number	Date of Draw	Bleed Number	RPR	VENDOR HIV Ag p24 S/Co	VENDOR ANTI HIV 1/2	VENDOR HBsAg	VENDOR ANTI HCV	TRINITY CAPTIA IgM S/Co	TRINITY CAPTIA IgG S/Co	Becton Dickinson RPR	FISHER Sure Vue RPR	INVERNESS FTA-ABS WAMPOLE SYPHILLIS
K-ZMC002-01	09/18/1996	N-10643-093	POS	NEG	NEG	NEG	NEG	0.29	2.49	Not Tested	POS	POS
K-ZMC002-02	07/11/1996	N-86246-062	POS	NEG	NEG	NEG	NEG	0.43	1.61	1:8	POS	POS
K-ZMC002-03	07/07/1996	N-75357-096	POS	NEG	NEG	NEG	NEG	0.50	1.35	1:2	POS	POS
K-ZMC002-04	07/08/1996	N-40538-063	POS	NEG	NEG	NEG	NEG	0.36	3.16	1:8	POS	POS
K-ZMC002-05	09/23/1996	N-56902-139	POS	NEG	NEG	NEG	NEG	0.36	2.95	Not Tested	POS	POS
K-ZMC002-06	09/24/1996	N-55123-020	POS	NEG	NEG	NEG	NEG	0.73	2.34	1:4	POS	POS
K-ZMC002-07	07/23/1996	G-54655-081	POS	NEG	NEG	NEG	NEG	1.02	3.28	1:8	POS	POS
K-ZMC002-08	08/19/1996	G-85547-014	POS	NEG	NEG	NEG	NEG	0.36	1.6	1:2	POS	INDETERM
K-ZMC002-09	10/13/1996	N-44833-014	POS	NEG	NEG	NEG	NEG	0.15	1.93	1:2	POS	POS
K-ZMC002-10	07/26/1996	N-13962-041	POS	NEG	NEG	NEG	NEG	1.24	3.7	Not Tested	POS	POS
K-ZMC002-11	06/28/1996	G-67031-036	POS	NEG	NEG	NEG	NEG	0.95	1.34	1:4	POS	POS
K-ZMC002-12	09/11/1996	N-45745-018	POS	NEG	NEG	NEG	NEG	1.17	2.27	Not Tested	POS	POS
K-ZMC002-13	10/09/1996	N-12969-079	POS	NEG	NEG	NEG	NEG	0.36	2.44	1:4	POS	POS
K-ZMC002-14	07/26/1996	N-92228-039	POS	NEG	NEG	NEG	NEG	0.95	4.54	Not Tested	POS	POS
K-ZMC002-15	08/07/1996	N-14210-108	POS	NEG	NEG	NEG	NEG	0.29	2.61	Not Tested	POS	POS