



GLOBAL HIV SURVEILLANCE PANEL
(Catalog No. K-GSP103)

LOT NUMBER: 0910-272-00069
For Research Use Only

MEMBER ID #	ORIGIN	DONOR REACTIVITY	ABBOTT HIV 1/2 Ab	CAMBRIDGE BIOTECH HIV-1 WESTERN BLOT	BIORAD MULTISPOT HIV 1	BIORAD MULTISPOT HIV 2	CHEMBIO HIV 1/2 STAT-PAK	CHEMBIO HIV 1/2 SURE CHECK	ORASURE ORAQUICK ADVANCE HIV 1/2	TRINITY UNI-GOLD RECOMBIGNEN HIV-1	MEDMIRA REVEAL G3 HIV-1
GSP103-01	USA	NORMAL	0.1	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG
GSP103-02	China	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-03	Ghana	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-04	Ghana	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-05	Thailand	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-06	USA	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-07	Cameroon	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-08	USA	NORMAL	0.1	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG
GSP103-09	Nigeria	HIV-1/2	>18	POS	POS	POS	POS	POS	POS	POS	POS
GSP103-10	China	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-11	USA	HIV-1/2	>18	POS	POS	POS	POS	POS	POS	POS	POS
GSP103-12	India	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-13	Ghana	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-14	Uganda	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-15	USA	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-16	USA	NORMAL	0.1	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG
GSP103-17	USA	HIV-1/2	>18	POS	POS	POS	POS	POS	POS	POS	POS
GSP103-18	Thailand	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-19	Ghana	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-20	India	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-21	Uganda	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-22	Cameroon	HIV-1	>18	POS	POS	NEG	POS	POS	POS	POS	POS
GSP103-23	Nigeria	HIV-1/2	>18	POS	POS	POS	POS	POS	POS	POS	POS
GSP103-24	USA	NORMAL	0.1	NEG	NEG	NEG	NEG	NEG	NEG	NEG	NEG

Run Date	16-Nov-09	03-Dec-09	24-Nov-09	24-Nov-09	24-Nov-09	24-Nov-09	24-Nov-09	24-Nov-09	24-Nov-09	24-Nov-09	24-Nov-09
Kit Lot #	79834M100	C0144-007	962473	962473	HIV090109	22090409	6617491	T224012	R3IA0118		
Exp. Date	27-Nov-09	24-May-09	29-Jun-10	29-Jun-10	12-Aug-11	20-Aug-11	20-Sep-10	14-Aug-10	31-Jul-10		

Data is presented as an S/Co ratio unless stated otherwise. Ratios >1.0 are considered to be reactive.

Notes on Data: This panel is for investigational use only and the data is offered for informational purposes only. Data is generated at ZeptoMetrix Corporation and by our Lab Partners in the United States and Europe. We will endeavor to update our data sheets whenever new test data becomes available to us. All specimens collected in this series are naturally occurring, unprocessed samples collected from HIV infected donors worldwide. No preservatives have been added.