

SPOTLIGHT

IMMUNO-TEK IgG/IgM Quantitative ELISAs




ZeptoMetrix Corporation
872 Main Street
Buffalo, NY 14202
Tel: 800 274 5487
Fax: 716 882 0959
Website:
www.ZeptoMetrix.com

IgG/IgM Quantitative ELISAs

Four assays for the measurement of human and mouse IgG and IgM antibodies are now available.

Each IMMUNO-TEK 96 test kit features a standard curve with a range from 7.8 to 125 ng/mL. The assays were designed for the measurement of antibodies in cell culture supernatants, ascites and other body fluids. These EIAs are especially useful in monitoring the production and purification of monoclonal antibodies.

The detector antibody in each kit has been specifically balanced to react uniformly with all subclasses of human or mouse IgG or IgM. All reagents necessary for performing each assay are included and provided in a liquid, ready-to-use format. or further product details.



IgG/IgM Quantitative ELISAs

The kits feature break-apart wells which allow for flexible batch sizes and helps minimize reagent waste. Total test time is less than 2 hours. The shelf life is 12 months from the date of manufacture.

Product inserts are available in PDF & HTML format on this web site. Please contact us directly for pricing and further details.

IMMUNO-TEK ELISAs

Human IgG Catalog [#08011](#)

Human IgM Catalog [#08011](#)



Mouse IgG Catalog [#0801180](#)

Mouse IgM Catalog [#0801181](#)

