



SPOTLIGHT

NOW AVAILABLE

**HIV-1 p24
Antigen ELISA
Bulk Kits**



ZeptoMetrix Corporation
872 Main Street
Buffalo, NY 14202
Tel: 800 274 5487
Fax: 716 882 0959
Website:
www.ZeptoMetrix.com

**HIV-1 p24
Antigen ELISA bulk kit**

ZeptoMetrix now has available a HIV-1 p24 Antigen ELISA bulk kit to meet the testing needs of large laboratories. The package design requires less storage space than single assays and provides a cost savings per kit for bulk purchases.

This assay is for *Research Use Only* and is intended for use in detection of HIV-1 p24 antigen in cell culture media as well as human sera and plasma. Each bulk kit contains five antibody coated microplates and all the necessary reagents to perform 480 determinations. The standard curve ranges from 7.8 pg/mL to 125 pg/mL.

The kit features break apart wells which allow for flexible batch sizes and minimal waste. Overnight and 2 hour protocols are provided.

Additional assays are available for immune complex dissociation and confirmation as well as extending the range of the standard curve for those samples with greater viral loads.

When ordering, please mention code WS-0402

