



Natural Human Interleukin-2 (IL-2) / T-Cell Growth Factor (TCGF)

Catalog #: 0801017

PRODUCT CHARACTERISTICS

Each lot of IL-2/TCGF is analyzed for its ability to stimulate proliferation of seven to ten day old PHA-transformed human T- lymphocytes. When used at a final concentration of approximately 10% (v/v), IL-2/TCGF will induce a minimum five-fold increase in cell concentration of such cells when seeded at a density of 2.0×10^5 / ml. IL-2/TCGF, a glycoprotein with a molecular weight of approximately 15,000 kDa, is produced from pooled human PHA-stimulated T-lymphocytes. IL-2/TCGF is purified by several chromatographic steps to remove PHA and interferon.

RECOMMENDED USAGE

This preparation will support growth of most IL-2/TCGF dependent lymphocyte cultures at a final concentration of 10% (v/v). The initial cell density should be approximately 2.0×10^5 / ml and the cell culture should be adjusted to this density every four to five days by additional feeding with 10% IL-2/TCGF. The initial cell population may need to be blast-transformed by either mitogenic or antigenic stimulation such as with PHA. Optimal concentration of IL-2/TCGF and/or mitogen must be determined experimentally by the investigator.

CONTENTS

Each 50 ml bottle contains approximately 25,000 BRMP (Biological Response Modifier Program) units of IL-2/TCGF at a concentration of approximately 500 BRMP units/ml. It is supplied as a sterile solution in 25mM HEPES buffered RPMI 1640 culture medium which is free of serum. IL-2/TCGF protein concentration ranges from 8-22 μ g/ml.



STORAGE AND STABILITY

The solution is stable at -10°C or below. The material may be re-frozen after thawing. Repetitive freezing and thawing is not recommended (aliquot if necessary). Thawed material may be stored at 4°C for short-term usage.

CAUTION: DO NOT USE IN HUMANS. FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The solution is stable at -10°C or below. The material may be re-frozen after thawing. Repetitive freezing and thawing is not recommended (aliquot if necessary). Thawed material may be stored at 4°C for short-term usage.

This product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.

REF	Catalog Number		Temperature Limitation
LOT	Lot Number	RUO	For Research Use Only
	Biological Risk		

ZeptoMetrix Corporation
878 Main Street
Buffalo, NY 14202
800-274-5487
www.zeptometrix.com