



***(BSL-3) St. Louis Encephalitis Virus
Strain: Parton
Infectious Culture Fluid (1 mL)**

Catalog Number: 0810080CF

PRODUCT DESCRIPTION:

St. Louis Encephalitis Virus (Strain: Parton) is an enveloped, icosahedral Flavivirus and contains a non-segmented, single-stranded positive-sense RNA.

Each frozen aliquot contains 1 mL of titered viral culture fluid.

INTENDED USE:

Viral culture fluids are sold as consumable testing materials; propagation or commercialization is prohibited without prior written consent from ZeptoMetrix. The suitability and performance characteristics should be determined by your laboratory for each intended usage.

These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The purchase of infectious microorganisms from ZeptoMetrix requires a Material Transfer Agreement (MTA).

BIOSAFETY:

St. Louis Encephalitis Virus is a Biosafety Level 3 microorganism and must be used within Biological Safety Level 3 facility. Please consult your institution's regulations regarding the use of this product. For a detailed discussion on biological safety see the 5th edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), published by the CDC at <http://www.cdc.gov/biosafety/publications/bmbl5/index.htm>.

PRECAUTIONS:

- Use Universal Precautions, this product is **potentially biohazardous**.
- Repetitive freezing and thawing is not recommended (aliquot material if necessary). Titer will be altered by multiple freeze-thaws.
- To avoid cross-contamination, use separate pipette tips for all reagents.

RECOMMENDED STORAGE:

Viral culture fluids should be stored at -65°C or below.

***EXPORT COMPLIANCE:**

This organism is currently included on the Commerce Control List (15 C.F.R. 774, Supplement 1). **Export or re-export of such organism, certain genetic elements thereof and/or genetically modified organism is prohibited without the appropriate authorization or license from the Bureau of Industry and Security per the U.S. Export Administration Regulations (15 C.F.R. 730-774).**

This product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.

REF	Catalog Number		Temperature Limitation
LOT	Lot Number		Expiration Date
RUO	For Research Use Only		Biological Risk
	Manufacturer		Date of Manufacture

ZeptoMetrix Corporation
878 Main Street
Buffalo, NY 14202
800-274-5487
www.zeptometrix.com