



***Human Metapneumovirus (hMPV-5) 5 Type B1
Strain: Peru3-2003
Infectious Culture Fluid (1 mL)**

Catalog Number: 0810158CF

PRODUCT DESCRIPTION:

Human Metapneumovirus (hMPV-5) 5 Type B1 (Strain: Peru3-2003) is an enveloped virus that contains a negative-sense single-stranded RNA.

Each frozen aliquot contains 1 mL of titered viral culture fluid.

INTENDED USE:

Viral culture fluids are sold as consumable testing materials; propagation or commercialization is prohibited without prior written consent from ZeptoMetrix. The suitability and performance characteristics should be determined by your laboratory for each intended usage.

These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The purchase of infectious microorganisms from ZeptoMetrix requires a Material Transfer Agreement (MTA).

BIOSAFETY:

Human Metapneumovirus is a Biosafety Level 2 microorganism and must be used within a Biological Safety Level 2 facility or cabinet. Please consult your institution's regulations regarding the use of this product. For a detailed discussion on biological safety see the 5th edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), published by the CDC at <http://www.cdc.gov/biosafety/publications/bmbl5/index.htm>.

PRECAUTIONS:

- Use Universal Precautions, this product is **potentially biohazardous**.
- Repetitive freezing and thawing is not recommended (aliquot material if necessary). Titer will be altered by multiple freeze-thaws.
- To avoid cross-contamination, use separate pipette tips for all reagents.

RECOMMENDED STORAGE:

Viral culture fluids should be stored at -65°C or below.

***Note: "This product is sold by ZeptoMetrix Corporation under license from ViroNovative B.V. under patent applications, including U.S. Patent Applications 10/371,099 and 10/371,122, and any patents that issue from applications related to PCT/NL02/00040 and PCT/US03/05271."**

This product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.

REF	Catalog Number	Temperature Limitation
LOT	Lot Number	Expiration Date
RUO	For Research Use Only	Biological Risk
Manufacturer		Date of Manufacture

ZeptoMetrix Corporation
878 Main Street
Buffalo, NY 14202
800-274-5487
www.zeptometrix.com