

**Varicella Zoster Virus  
Isolate: A  
Infectious Culture Fluid (1 mL)**

## PRODUCT DESCRIPTION:

Varicella Zoster Virus (VZV) (Isolate: A) is a Human Herpes Virus (HHV-3) with a diameter of approximately 200 nm. It is an enveloped icosahedral virus that contains a nucleocapsid and double-stranded linear DNA.

Each frozen aliquot contains 1 mL of viral culture fluid.

Varicella Zoster Virus does not always demonstrate distinct CPE in culture. Therefore, VZV culture fluid is assayed by PCR (copies/mL).

## INTENDED USE:

Viral culture fluids are sold as consumable testing materials; propagation or commercialization is prohibited without prior written consent from ZeptoMetrix. The suitability and performance characteristics should be determined by your laboratory for each intended usage.

These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

**The purchase of infectious microorganisms from ZeptoMetrix requires a Material Transfer Agreement (MTA).**

## BIOSAFETY:

Varicella Zoster Virus is a Biosafety Level 2 (BSL-2) microorganism and must be used within a BSL-2 facility in a biosafety cabinet (BSC). Please consult your institution's regulations regarding the use of this product. For a detailed discussion on biological safety see the 5<sup>th</sup> edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), published by the CDC at <http://www.cdc.gov/biosafety/publications/bmbl5/index.htm>.








## PRECAUTIONS:

- Use Universal Precautions, this product is **potentially biohazardous**.
- Repetitive freezing and thawing is not recommended (aliquot material if necessary). Titer may be altered by multiple freeze-thaws.
- To avoid cross-contamination, use separate pipette tips for all reagents.

## RECOMMENDED STORAGE:

Viral culture fluids should be stored at -65°C or below.

PI0810172CF  
Revision: 05  
Effective Date: 08/05/2021

 REF	Catalog Number	 TEMP	Temperature Limitation
 LOT	Batch Code	 EXPI	Expiration Date
 RUO	For Research Use Only	 BIO	Biological Risk
 MFG	Manufacturer		

[www.ZeptoMetrix.com](http://www.ZeptoMetrix.com)

ZeptoMetrix LLC • 878 Main Street, Buffalo, NY 14202 USA • Tel (800) 274-5487

This product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.