



HCV Seroconversion Panel Donor No. 22708D
Product Code: HCV9094

Donor Demographics: Donor source plasma was collected from a U.S. licensed blood collection facility.
Donor was un-treated and asymptomatic throughout the collection period.



Panel Member	Donor Number	HCV	QUEST HCV Ab	ROCHE HCV Ampliscreen	ROCHE HCV Monitor
9094-01	22708D-01	Negative	0.28	Negative	Not Detected
9094-02	22708D-02	Negative	0.34	Negative	Not Detected
9094-03	22708D-03	Negative	1.86	Negative	Not Detected
9094-04	22708D-04	Positive	2.32	Negative	Not Detected
9094-05	22708D-05	Positive	3.00	Negative	Not Detected

Notes on Data:

1. This panel is for **Research Use Only. Not for use in diagnostic procedures.**
2. Data is generated at ZeptoMetrix™ Corporation and by our Laboratory Partners in the United States and Europe.
3. All specimens collected in this longitudinal series are unadulterated 4% sodium citrated plasma samples collected from a single donor in the United States.
4. No preservatives have been added.

This Product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.

ZeptoMetrix Corporation • 25 Kenwood Circle, Franklin, MA 02038 • Tel (508) 553-5800 • Fax (508) 520-1525
For Customer Support, please visit www.zeptometrix.com or email custserv@zeptometrix.com