



**HCV Seroconversion Panel Donor No.1038851**  
**Product Code: HCV9095**

**Donor Demographics:** Donor source plasma was collected from a U.S. licensed blood collection facility.  
Donor was un-treated and asymptomatic throughout the collection period.



Panel Member	Donor Number	HCV	DIASORIN HBsAb	QUEST HCV Ab	ROCHE HCV Ampliscreen	ROCHE HCV Monitor
9095-01	1038851-01	Negative	Not Tested	0.04	<b>Positive</b>	<b>159,000</b>
9095-02	1038851-02	Negative	Not Tested	0.04	<b>Positive</b>	<b>54,000</b>
9095-03	1038851-03	Negative	Not Tested	<b>11.70</b>	<b>Positive</b>	<b>141,000</b>
9095-04	1038851-04	<b>Positive</b>	Not Tested	<b>18.10</b>	<b>Positive</b>	<b>17,500</b>
9095-05	1038851-05	<b>Positive</b>	0.08	<b>18.40</b>	<b>Positive</b>	<b>42,600</b>

**Notes on Data:**

1. This panel is for **Research Use Only. Not for use in diagnostic procedures.**
2. Data is generated at ZeptoMetrix™ Corporation and by our Laboratory Partners in the United States and Europe.
3. All specimens collected in this longitudinal series are unadulterated 4% sodium citrated plasma samples collected from a single donor in the United States.
4. No preservatives have been added.

**This Product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.**

ZeptoMetrix Corporation • 25 Kenwood Circle, Franklin, MA 02038 • Tel (508) 553-5800 • Fax (508) 520-1525  
For Customer Support, please visit [www.zeptometrix.com](http://www.zeptometrix.com) or email [custserv@zeptometrix.com](mailto:custserv@zeptometrix.com)