

# SAFETY DATA SHEET (SDS)

## 1. Product Identification:

<b>Product Name:</b>	NATtrol™ SARS-Related Coronavirus 2 Positive Control
<b>Catalog Number:</b>	NATSARS(COV2)-ERC-IVD
<b>Product ID Number (CAS #)</b>	No Data Available
<b>Intended use of the product:</b>	Class 1 exempt, unassayed, <i>in vitro</i> , diagnostic external run controls intended for use with qualitative molecular assays
<b><u>ZeptoMetrix LLC</u></b>	
25 Kenwood Circle	
Franklin MA 02038	
Phone: (508)553-5800	
Fax: (508)553-5899	
<b>Emergency telephone number:</b>	1-800-274-5487 (US Toll Free) Monday-Friday 8:15am – 4:30pm EST

## 2. Hazard(s) Identification:

<b>Classification of the substance or chemical:</b>	No Data Available
<b>Classification:</b>	Not a hazardous substance or mixture
<b>Hazard Pictograms:</b>	None
<b>Signal word:</b>	Warning
<b>Hazard statements:</b>	Potential Biohazard, Use Universal Precautions
<b>Precautionary statement:</b>	May present unknown hazards and should be handled as if capable of transmitting infectious agents.
<b>Other hazards:</b>	Contains material derived from human sources and may be considered a biohazard and/or Regulated Medical Waste in your state. Check your Local Environmental Regulations.
<b>Unknown acute toxicity:</b>	No Data Available

## 3. Composition/Information on Ingredients:

<u>Component Name</u>	<u>CAS Number</u>	<u>Concentration</u>	<u>Classification</u>
Trade Secret- Undisclosed Identity	Not Applicable	Not Applicable	Not Applicable

## 4. First Aid/Emergency Measures:

<b><u>Description of first aid measures</u></b>	
First-aid measures general:	If you feel unwell, seek medical advice.
First-aid measures after Eye Contact:	Rinse cautiously with copious amounts of water for several minutes.
First-aid measures after Skin Contact:	Remove contaminated clothing and immediately rinse skin with copious amounts of water followed by washing with soap and copious amounts of water.
First-aid measures after Ingestion:	If the person is unconscious, seek emergency medical attention: never give anything by mouth to an unconscious person. If the person is conscious, wash mouth out with copious amounts of water and call a physician. Do not induce vomiting unless directed to do so by a physician.
First-aid measures after Inhalation:	If person is unconscious, seek emergency medical attention. If person is conscious, remove to fresh air and call a physician.
<b><u>Most important symptoms and effects, both acute and delayed</u></b>	
Symptoms/injuries after Eye Contact:	Seek medical assistance if required.
Symptoms/injuries after Skin Contact:	Seek medical assistance if required.
Symptoms/injuries after Ingestion:	Seek medical assistance if required.
Symptoms/injuries after Inhalation:	Seek medical assistance if required.
<b><u>Indication of any immediate medical attention and special treatment needed:</u></b> No data available	

## 5. Fire and Explosion Measures:

<b>Extinguishing Media</b>	
Suitable extinguishing media:	Use an extinguishing agent suitable for the surrounding area.
Unsuitable extinguishing media:	Use an extinguishing agent suitable for the surrounding area.
<b>Special hazards arising from the substance or mixture</b>	
Fire hazard:	No data available
Explosion hazard:	No data available
Reactivity:	No data available
<b>Advice for firefighters</b>	
Firefighting instructions:	Exercise caution when fighting a fire.
Protection during firefighting:	Firefighters should wear protective gear. Do not enter fire area without proper protective equipment, including respiratory protection.

## 6. Accidental Release Measures:

<b>Personal precautions, protective equipment and emergency procedures</b>	
General measures:	Use appropriate personal protective equipment (PPE) and appropriate laboratory procedures .
<b>For non-emergency personnel</b>	
Protective equipment:	Use appropriate personal protective equipment (PPE).
Emergency procedures:	Use appropriate emergency procedures.
<b>For emergency responders</b>	
Protective equipment:	Use appropriate personal protective equipment (PPE).
Emergency procedures:	Use appropriate emergency procedures.
<b>Methods and material for containment and cleaning up</b>	
For containment:	Contain material in accordance to State and Federal regulations.
Methods for cleaning up:	Follow and dispose of material in accordance to State and Federal waste disposal regulations.

## 7. Handling and Storage:

<b>Precautions for safe handling</b>	
Precautions for safe handling:	Handle in accordance with good laboratory practices and Biosafety Level 2 safety procedures. Use Universal Precautions.
<b>Conditions for safe storage, including any incompatibilities</b>	
Storage conditions:	Store in the original container. Keep container closed when not in use.
Incompatible products:	No data available
Incompatible materials:	No data available
Storage:	Recommended storage temperature 2-8°C (36-46°F).

## 8. Exposure Controls and Personal Protection:

<b>Exposure Controls</b>	
Appropriate engineering controls:	Handle in accordance with good laboratory practices and safety procedures.
Personal Protective Equipment:	Use laboratory coat, protective gloves, safety glasses and suitable protective clothing.
Pictograms:	
Eye/Face protection:	Wear in accordance with good laboratory practices and safety procedures.
Skin and body protection:	Wear in accordance with good laboratory practices and safety procedures.
Respiratory protection:	Wear in accordance with good laboratory practices and safety procedures.
Hand protection:	Wear in accordance with good laboratory practices and safety procedures.
Mechanical protection:	Wear in accordance with good laboratory practices and safety procedures.
Special work practices:	Wear in accordance with good laboratory practices and safety procedures.

## 9. Physical and Chemical Properties:

<b>Information on basic physical and chemical properties</b>	
Physical State:	Liquid
Appearance/color:	Clear, pale amber solution
Odor:	No data available
Odor threshold:	No data available
pH:	No data available
Melting point/freezing point:	No data available
Boiling point:	No data available
Flash point:	No data available
Evaporation rate:	No data available
Flammability:	No data available
Upper/lower flammability or explosive limits:	No data available
Vapor pressure:	No data available
Vapor density:	No data available
Relative density:	No data available
Solubility(ies):	No data available
Partition coefficient:	No data available
Auto-ignition temperature:	No data available
Decomposition temperature:	No data available
Viscosity:	No data available
<b>Other information:</b>	No data available

## 10. Stability and Reactivity:

Reactivity:	No data available
Chemical stability:	Stable under recommended storage conditions.
Possibility of hazardous reactions:	No data available
Conditions to avoid:	No data available
Incompatible materials:	No data available
Hazardous decomposition products:	No data available

## 11. Toxicological Information:

<b>Information on toxicological effects</b>	
Skin corrosion/irritation:	May cause mild irritation.
Eye damage/irritation:	May cause mild irritation.
Respiratory damage/irritation:	May cause mild irritation.
Ingestion damage/irritation:	May cause mild irritation.
Specific target organ toxicity (single exposure):	No data available
Specific target organ toxicity (repeated exposure):	No data available
Numerical measure of toxicity:	No data available
Symptoms/injuries after skin contact:	No data available
Symptoms/injuries after eye contact:	No data available
Symptoms/injuries after inhalation:	No data available
Symptoms/injuries after ingestion:	No data available

## 12. Ecological Information:

Ecotoxicity:	The product does not contain significant quantities of ingredients that are environmentally toxic.
Persistence and degradability:	The product is unlikely to persist in the environment.
Bioaccumulative potential:	No data available
Mobility in soil:	No data available
Other adverse effects:	No data available

## 13. Disposal Information:

Waste disposal recommendations:	Dispose of material in accordance to local, State and Federal waste disposal regulations concerning potentially biohazardous waste.
---------------------------------	---

## 14. Transport Information:

In accordance with IATA Regulations	
UN number	No data available
UN proper shipping name	No data available
Department of Transportation (DOT)	No data available
Hazard Classes	No data available
Hazard labels (DOT)	No data available
Packing group (DOT)	No data available
Environmental hazards	No data available
Transport in bulk	No data available
Special precautions	No data available

## 15. Regulatory Information:

SARA 311/312: Hazard categories for SARA Section 311/312 Reporting	No data available
Canadian WHMIS Classification:	No data available
EU Classification (90/492/EE):	No data available
EU Risk and Safety Phrases:	No data available
California Proposition 65:	This product does not contain chemicals known to the State of California to cause cancer, birth defects or any other reproductive harm.

## 16. Other Information:

Preparation Date	07/08/2022
------------------	------------

### **DISCLAIMER**

The information, data and recommendations contained herein ("Information") are based upon information believed by ZeptoMetrix LLC ("ZM") to be reliable, however ZM neither warrants the accuracy or completeness of the Information. ZM MAKES NO REPRESENTATION OR WARRANTY, EITHER EXPRESS OR IMPLIED, WHATSOEVER, WITH REGARD TO THE INFORMATION OR THE PRODUCT DESCRIBED HEREIN, INCLUDING ANY WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR USE, NON-INFRINGEMENT. , THAT THE PRODUCT WILL NOT POSE A SAFETY OR HEALTH RISK, OR THAT THE PRODUCT WILL ACCOMPLISH ANY PARTICULAR RESULT. ZM shall not be held liable or responsible for any damages or losses resulting from use of the product described herein or the Information.

All materials and mixtures may present unknown hazards and should be used with caution. When necessary or appropriate, independent opinion regarding the risk of handling or exposure should be obtained from trained professionals.